


TOPTICA's role in ATI Q

The ATI Q project will produce three quantum computing demonstrators based on different ion trap technologies. For each demonstrator, laser light with well-controlled frequency, intensity, and switching capabilities is needed to manipulate the experimental system and perform the essential operations of a quantum computer. TOPTICA will leverage a modular platform designed for light processing to extend this platform to include the wavelengths and functionalities necessary to produce rack-mounted light processing units adapted to the needs of each demonstrator. In addition, TOPTICA will contribute a rack-mounted solution that prepares laser light for qubit addressing of multiple ion zones for one of the demonstrators.

Laser Beam Conditioning




Full solution for laser beam conditioning in a rack mounted solution, including:

- Beam routing and switching
- Frequency shifting, sidebands (AOM/EOMs)
- Phase and polarization control

Modular platform: Light Processing Units

Prototype designed in project IQuAn (quantum processor):



Light Processing Unit (LPU) – example configuration

Modules

I. Interface Cube

- Polarization cleaning
- Beam splitting, folding
- Power monitoring


II. Double pass AOM

- Frequency shifts
- Intensity modulation
- Switching

III. Switch EOM

IV. Sideband EOM

V. Automated Fiber Coupling




Rack based solution:

- Compact and hands-off
- Fully fiber coupled
- Modular construction
- Passive stability
- Cable management
- Thermal management

Expansion of wavelengths and functionalities

Goal of research within ATI Q project:

- Extension beyond Ca⁺ wavelengths
- Investigation of added functionalities
 - Shutters
 - Beam combining
- Diverse configurations
 - 3 demonstrator approaches
 - 4 species
- Delivery of **25-30 custom LPUs** in total



Ion	Wavelengths in use
Yb ⁺	369 nm, 399nm, 760 nm, 935 nm
Ba ⁺	413nm, 493nm, 614 nm, 650 nm, 1762 nm
Sr ⁺	405 nm, 422 nm, 461 nm, 1033 nm, 1091 nm
Ca ⁺	375 nm, 397 nm, 422 nm, 729 nm, 854 nm, 866 nm



Ca⁺, Sr⁺




Ca⁺, Sr⁺



Yb⁺, Ba⁺

Multiple-zone Raman addressing unit

Goal: Manipulation of ground state qubits with Raman transitions for quantum processor in Mainz (JGUM)



- Up to 16 addressing beams
- Frequency, intensity, polarization, phase control
- High passive stability required for performing gates



SPONSORED BY THE

